

Pin Covers

Our Pin Covers help to protect patients, staff, and clothing from the sharp ends of Steinmann Pins and K-Wires. These are made of thermoplastic and are available in three color-coded sizes to fit all eleven sizes of Steinmann Pins and K-Wires. They are latex-free, sterile, and packaged for single-use.

- Yellow (Small) – Packaged 1/Pouch to match applications that involve one exposed end.
- Blue (Medium) and Red (Large) – Packaged 2/Pouch to match applications that involve two exposed ends.



Part No.	Color	Fits This Size Range/Inches	Quantity Per Box
101001pbx	Yellow	0.028 – 0.062	1/ Pouch, 60 Pouches/Box
101002pbx	Blue	0.078 – 0.125	2/Pouch, 30 Pouches/Box
101003pbx	Red	0.141 – 0.188	2/Pouch, 30 Pouches/Box
101000pbx	Assorted	Assorted	50 Assorted Pouches/Box

Disposable Light Handles

We offer two systems for convenience: Rigid Light Handle or Flexible Cover. Eliminate costly reprocessing with these easy to use disposable Rigid Handles or Flexible Covers. All are sterile, latex-free, and intended for single-use.

Adaptors are available for most of the popular lights. Please call for further adapter information.

Rigid Handle




Figure A:
Illustrates our Rigid Light Handle. Our threaded design fits most adaptors made.

Flexible Cover




Figure B:
The Flexible Cover was designed using our unique push system allowing you to easily apply the cover. The pull tab makes for fast removal.

Part No.	Description	Quantity Per Box
LT-R01	Rigid Handle	24
LT-F01b	Flexible Cover (1 Each)	80
LT-F02	Flexible Cover (2 Each)	80
LT-F03	Flexible Cover (3 Each)	80