

### Instrument Identification Tape

Tape makes organizing and tracking valuable instruments easy! Tape is flexible, adheres securely to instruments, and will not chip, flake, peel, or crack — and with 75+ colors, symbols, and character options, the color-coding options are nearly endless. If a change is needed, simply remove the ColorTrack I.D. Tape and replace with your new color choice. Latex-free.

#### Sheet Tape

Each 8½" x 11" sheet has 374 total inches of tape and contains:

- 44 each of 2 inch strips
- 44 each of 1½ inch strips
- 220 each of 1 inch strips

#### Solid Color

¼" Height / Various widths

Part No.	Description
152001eea	White
152002eea	Royal Blue
152003eea	Kelly Green
152004eea	Red
152005eea	Yellow
152006eea	Purple
152007eea	Orange
152008eea	Brown
152009eea	Black
152010eea	Gray
152011eea	Lavender
152012eea	Beige
152013eea	Light Blue
152014eea	Lime Green
152015eea	Pink
152017eea	Plum
152018eea	Mint Green

#### Dots

¼" Height / Various widths

Part No.	Description
158040eea	Dots, Fuchsia
158041eea	Dots, Medical Blue
158042eea	Dots, Emerald Green
158043eea	Dots, Violet

#### Stars

¼" Height / Various widths

Part No.	Description
158030eea	Stars, Fuchsia
158031eea	Stars, Medical Blue
158032eea	Stars, Emerald Green
158033eea	Stars, Violet

#### White Stripe

¼" Height / Various widths

Part No.	Description
154002eea	White/Royal Blue
154003eea	White/Kelly Green
154004eea	White/Red
154005eea	White/Yellow
154007eea	White/Orange
154008eea	White/Brown
154011eea	White/Lavender
154013eea	White/Light Blue
154014eea	White/Lime Green
154015eea	White/Pink
154017eea	White/Plum

#### Diagonal Stripes

¼" Height / Various widths

Part No.	Description
159001eea	Diagonal, Purple/Black
159002eea	Diagonal, Orange/Black
159003eea	Diagonal, Blue/Black
159004eea	Diagonal, Rubine/Black
159005eea	Diagonal, Yellow/Black
159006eea	Diagonal, Green/Black
159007eea	Diagonal, White/Black

#### Tri-Color Stripe

⅜" Height / Various widths

Part No.	Description
155051eea	Red/White/Blue
155052eea	Red/White/Black
155053eea	Red/Yellow/Black
155054eea	Red/Black/Green

#### Black Stripe

¼" Height / Various widths

Part No.	Description
153001eea	Black/White
153003eea	Black/Kelly Green
153004eea	Black/Red
153005eea	Black/Yellow
153007eea	Black/Orange
153010eea	Black/Gray
153011eea	Black/Lavender
153012eea	Black/Beige
153013eea	Black/Light Blue
153014eea	Black/Lime Green
153015eea	Black/Pink
153017eea	Black/Plum

#### Triangles

¼" Height / Various widths

Part No.	Description
158050eea	Triangles, Fuchsia
158051eea	Triangles, Medical Blue
158052eea	Triangles, Emerald Green
158053eea	Triangles, Violet

#### Numbers

¼" Height / Various widths

Each sheet contains numbers 0–8

Part No.	Description
11111	158020eea Numbers, Fuchsia
22222	158021eea Numbers, Medical Blue
33333	158022eea Numbers, Emerald Green
44444	158023eea Numbers, Violet

**Assorted Symbols**

¼" Height / Various widths

Part No.	Description
158060eea	Red Hearts on White
158070eea	Smiley Face, Yellow/Black
158080eea	Alphabet, Fuschia/Wt, A-E
158061eea	White Hearts on Red

**4-Color**

¼" Height / Various widths

Part No.	Description
152031eea	Red/Yellow/Purple/Orange
152032eea	Blue/White/Black/Green

**Specialty Stripe**

¼" Height / Various widths

Part No.	Description
156003eea	Blue/Green
156004eea	Blue/Red
157002eea	Yellow/Blue
157003eea	Yellow/Green
157004eea	Yellow/Red

**Code to Specialty**

¼" Height / Various widths

Part No.	Description
ENT 158090eea	Mint/White ENT
ER ER 158091eea	Purple/White ER
GEN 158092eea	Royal Blue/White GEN
OB/GYN 158093eea	Lt Blue/White OB/GYN
ORTHO 158094eea	Kelly Green/White ORTHO
NEURO 158095eea	Orange/White NEURO
CARDIAC 158096eea	Red/White CARDIAC

**Roll Tape**

Roll Tape ¼" x 300"

Part No.	Description
151001eea	White
151002eea	Royal Blue
151003eea	Kelly Green
151004eea	Red
151005eea	Yellow
151006eea	Purple
151007eea	Orange
151008eea	Brown
151009eea	Black

Roll Tape ⅛" x 250"

Part No.	Description
151011eea	White
151012eea	Royal Blue
151013eea	Kelly Green
151014eea	Red
151015eea	Yellow
151016eea	Purple
151017eea	Orange
151018eea	Brown
151019eea	Black



**Application Instructions:**

For secure placement through multiple processing cycles, ColorTrack I.D. Tape should be wrapped 1½ times around the Instrument and excess length should be trimmed. Instrument sterility under ColorTrack I.D. Tape is assured after processing with properly applied tape.\*

\* Test Results Available On Request

Size and color of tapes shown may vary from actual products.

**Same Day Tape Remover**

An easy, instrument friendly solution to remove ColorTrack Tape. Eliminates the risk of cuts and sticks common with scalpel or utility knife tape removal. Neutral scent makes Same Day Tape Remover comfortable to use.

Part No.	Description
151066eea	1 Gallon Container

