

**Surgical Skin Markers**

Our surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery. Pens contain prep-resistant gentian violet ink and are available with or without ruler and blank labels. The only exception to this is our Taper Tip All-Skin Marker which contains a deeper blue ink for demarcation on darker skin types. A variety of barrel and tip options are available depending on your need and preference. Single-use and latex-free.

**Classic Fine Tip Skin Marker**



Part No.	Description	Qty
2530	Marker, Ruler, and Label Set — Sterile	12/Box
2531	Marker — Sterile	12/Box
2630	Marker, Ruler, and Label Set — Sterile	50/Box
2631	Marker — Sterile	50/Box

**Classic Regular Tip Skin Marker**



Part No.	Description	Qty
2550	Marker, Ruler, and Label Set — Sterile	12/Box
2551	Marker — Sterile	12/Box
2650	Marker, Ruler, and Label Set — Sterile	50/Box
2651	Marker — Sterile	50/Box

**Modern Regular Tip Skin Marker**



Part No.	Description	Qty
2712	Marker, Ruler, and Label Set — Sterile	12/Box
2750	Marker, Ruler, and Label Set — Sterile	50/Box
2751	Marker — Sterile	50/Box

**Taper Tip Skin Marker**



Part No.	Description	Qty
2730	Marker, Ruler, and Label Set — Sterile	50/Box
2731	Marker — Sterile	50/Box

**Taper Tip All-Skin Marker**



Part No.	Description	Qty
2756	Marker, Ruler, and Label Set — Sterile	50/Box
2756BN	Marker — Non-Sterile	500/Case

**Securiline® Taper Tip Skin Marker**



Part No.	Description	Qty
1000-00-PDG	Marker and Ruler Set — Sterile	10/Box, 10 Boxes/Case
1002-00-PDG	Marker and Ruler Set — Sterile	100/Box
100L-00-PDG	Marker, Ruler, and Label Set — Sterile	100/Box
1003-00-PDG	Marker and Ruler Set — Non-Sterile	100/Box

**WRITESITE® Preoperative Marking Products**

Our non-sterile markers or stamp pad with gentian violet ink are designed to be used in your preoperative patient surgical site marking protocol. The newest addition to this line, the WRITESITE® Jr., generates less waste for this single-use item. Both markers are offered in our regular tip. Single-use and latex-free.

**WRITESITE® Preoperative Surgical Site Marker**



Part No.	Description	Qty
2753BN	Marker — Non-Sterile	100/Box

**WRITESITE® Jr. Small Preoperative Skin Marker**



Part No.	Description	Qty
2900BN	Marker — Non-Sterile	100/Box

**WRITESITE® Preoperative Surgical Site "YES" Stamp**



Part No.	Description	Qty
WS-001NS	"Yes" Stamp — Non-Sterile	75/Box

**Utility Markers**

Our utility markers provide a permanent method of identifying instruments, specimens, and other objects. The markers have black ink, are waterproof, non-smearing, and can be used on most surfaces. Utility markers are available with or without nine black and white or twenty color medication labels and come in regular or taper tip depending on your need. We also offer the Multi-Ink marker, which has gentian violet skin ink on one end and black utility ink on the other. Latex-free.

**Regular Tip Utility Marker**



Part No.	Description	Qty
2850	Marker and 9 Label Set — Sterile	50/Box
2850BN	Marker and 9 Label Set — Non-Sterile	100/Box
2851	Marker — Sterile	50/Box
2851BN	Marker — Non-Sterile	100/Box

**Taper Tip Utility Marker**



Part No.	Description	Qty
2830	Marker and 9 Label Set — Sterile	50/Box
2830BN	Marker and 9 Label Set — Non-Sterile	100/Box
2831	Marker — Sterile	50/Box
2831BN	Marker — Non-Sterile	100/Box
ML5832	Marker and 20 Medication Label Set — Sterile	50/Box

**Multi-Ink Marker**



Part No.	Description	Qty
2754	Marker and 9 Label Set — Sterile	12/Box
2754BN	Marker and 9 Label Set — Non-Sterile	100/Box
2755	Marker — Sterile	12/Box
2755BN	Marker — Non-Sterile	100/Box

Select markers are available with Blank Labels, Rulers, or Medication Labels as noted.



Blank Labels



Rulers



20 Color Medication Labels



9 Black and White Medication Labels

**Securline® Lab Markers**

Our lab markers are versatile and can be used to make permanent markings on plastic, metal, and porcelain labware. The unique ink has the ability to write on labware even when the surface is damp. These markers are available in both fine and extra-fine tip as well as in black, blue, green, or red ink. Latex-free.

**Securline® Extra-Fine Tip Lab Marker**



Part No.	Description	Qty
1300-13-PDC	Blue Marker — Non-Sterile	10/Box
1300-15-PDC	Green Marker — Non-Sterile	10/Box
1300-16-PDC	Red Marker — Non-Sterile	10/Box
1300-20-PDC	Black Marker — Non-Sterile	10/Box

**Securline® Fine Tip Lab Marker**



Part No.	Description	Qty
1400-13-PDC	Blue Marker — Non-Sterile	10/Box
1400-15-PDC	Green Marker — Non-Sterile	10/Box
1400-16-PDC	Red Marker — Non-Sterile	10/Box
1400-20-PDC	Black Marker — Non-Sterile	10/Box

**Securline® Superfrost Marker**

The Superfrost Marker has permanent ink that resists common solvents used in histology and cytology. This pen has an extra-fine tip and is available in black or red ink. Latex-free.

**Securline® Extra-Fine Tip Superfrost Marker**



Part No.	Description	Qty
1450-16-PDC	Red Marker — Non-Sterile	10/Box
1450-20-PDC	Black Marker — Non-Sterile	10/Box