

**Bulldog Vascular Clamps**

Aspen surgical provides three varieties of Bulldog Vascular Clamps for use in occluding vessels, depending on desired clamping force, jaw opening, and jaw length. All are sterile, single-use, and latex-free.

**Atraumatic Bulldog Clamp**

The jaws of the Atraumatic Bulldog Clamp are coated with a soft material that gently grips the vessel reducing trauma to the patient, while still successfully occluding the vessel. The clamps have virtually the same clamping force of competitors while requiring less force to apply.

- Clamping Force – 1.048 lbs
- Jaw Opening – 6mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031001pbx	Atraumatic Bulldog Clamp	5 Pair

**Secure-Hold Bulldog Clamp**

The Secure-Hold Bulldog Clamp utilizes a soft material on the jaws of the clamp to occlude the vessel during use. This clamp is designed much like the Atraumatic Clamp to reduce a patient's trauma, while still occluding the vessel.

- Clamping Force – 1.149 lbs
- Jaw Opening – 6mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031003pbx	Secure-Hold Bulldog Clamp	5 Pair

**Extra-Traction Bulldog Clamp**

The Extra-Traction Bulldog Clamp has a serrated design to occlude larger and tougher vessels. The clamp has a larger jaw opening than the competitor, while still maintaining similar opening and clamping forces.

- Clamping Force – 1.556 lbs
- Jaw Opening – 10mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031004PBX	Extra-Traction Bulldog Clamp	5 Pair

**Syringe and Scalpel Handling System**

The Syringe and Scalpel Handling System is used to provide safe and easy handling and storage for syringes and scalpels. This item consists of an adhesive-backed block that firmly grips and holds the syringe cap, syringe, or scalpel. It is available with slots for one, two, or three syringes or scalpels, is intended for single use, and is latex-free.



Part No.	Block Size	Width mm	Width Inches	Quantity Per Box
<b>Sterile</b>				
111001pbx	Single	31.8	1.25	5
111002pbx	Double	44.5	1.75	5
111003pbx	Triple	69.9	2.75	5
<b>Non-Sterile</b>				
111001bbg	Single	31.8	1.25	12
111002bbg	Double	44.5	1.75	12
111003bbg	Triple	69.9	2.75	12