

**SafeClean™ Instrument and Channel Cleaning Brushes** are designed to clean surgical/medical instruments and endoscopic devices. The complete line of SafeClean™ Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices. Latex-free.

#### General Purpose Brushes

General Purpose Brushes can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241007bbg	1mm (0.039")	18"	0.5"	2
241010bbg	1.6mm (0.063")	18"	0.5"	2
241015bbg	1.8mm (0.071")	18"	0.5"	2
241019bbg	2mm (0.079")	12"	1"	2
241020bbg	2.2mm (0.086")	18"	0.5"	2
241024bbg	2.5mm (0.098")	24"	2.5"	2
241025bbg	2.7mm (0.106")	18"	1"	2
241029bbg	3mm (0.118")	12"	1"	2
241030bbg	3.2mm (0.126")	18"	2"	2
241032bbg	4mm (0.158")	12"	1"	2
241033bbg	5mm (0.197")	12"	2.5"	2
241034bbg	5mm (0.197")	24"	2.5"	2
241035bbg	5mm (0.197")	16"	2.5"	2
241040bbg	6mm (0.236")	16"	2.5"	2
241041bbg	6mm (0.236")	24"	2.5"	2
241042bbg	7mm (0.275")	16"	2.5"	2
241043bbg	7mm (0.275")	24"	2.5"	2
241044bbg	10mm (0.393")	16"	2.5"	2
241045bbg	11mm (0.433")	16"	2.5"	2
241046bbg	11mm (0.433")	24"	2.5"	2
241047bbg	10mm (0.393")	24"	2.5"	2
241050bbg	12mm (0.472")	16"	2.5"	2
241055bbg	15mm (0.590")	16"	2.5"	2
241060bbg	20mm (0.787")	16"	2.5"	2

#### Kits

Make ordering easy with all of the brushes you need in one kit!

##### Small Cleaning Brush Kit

This kit includes various small general purpose cleaning brushes. Latex-free.

##### Part No.

**241100bbg — 1-Each**

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

##### Large Cleaning Brush Kit

This kit includes various large general purpose cleaning brushes. Latex-free.

##### Part No.

**241200bbg — 1-Each**

241035bbg

241040bbg

241045bbg

241050bbg

241055bbg

241060bbg

##### "Best Sellers" Cleaning Brush Kit

This kit includes a selection of our most popular general purpose cleaning brushes. Latex-free.

##### Part No.

**241300bbg — 1-Each**

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

241035bbg

241040bbg

241045bbg

241050bbg

241055bbg

241060bbg

##### General Cleaning Brush Kit

This kit contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes. Latex-free.

##### Part No.

**241600bbg — 3-Each**

241025bbg

241030bbg

241035bbg

241040bbg

241041bbg

241045bbg

241046bbg

241050bbg

241055bbg

241060bbg

242001bbg

242002bbg

242003bbg

243001bbg

243004bbg

##### Micro/Ophthalmic Brush Kit

This kit contains a selection of our smallest general purpose cleaning brushes. Latex-free.

##### Part No.

**241700bbg — 3-Each**

241007bbg

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

**Suction Tube Cleaning Brushes**

Suction Tube Cleaning Brushes are designed to fit into and effectively clean the smallest bore devices and instruments used in surgery, MIS, and therapy procedures today. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241003bbg	3 French (1mm, 0.039")	6"	0.5"	2
241005bbg	5 French (1.67mm, 0.066")	6"	0.5"	2
241006bbg	6 French (2mm, 0.079")	8"	0.5"	2
241008bbg	7 French (2.3mm, 0.092")	8"	0.5"	2
241009bbg	8 French (2.7mm, 0.105")	8"	0.5"	2
241011bbg	9 French (3mm, 0.118")	8"	0.5"	2
241012bbg	10 French (3.3mm, 0.131")	8"	0.5"	2
241013bbg	11 French (3.7mm, 0.144")	8"	0.5"	2
241014bbg	12 French (4mm, 0.158")	8"	0.5"	2

**Fan-Tip Brush**

Fan-Tip Brushes feature bristles that wrap around the end of the brush for superior cleaning of bottles, jars, and other large instruments. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
242001bbg	4.8mm-12.7mm Tapered (0.189"-0.5")	7"	3"	2
242002bbg	13mm (0.512")	6"	3"	2
242003bbg	16mm (0.623")	5"	2.5"	2
242046bbg	9.6mm (0.384")	24"	2.5"	2

**Tracheostomy Brushes**

Use the right brush for the right job. Our Tracheostomy Brushes are uniquely designed to clean tubes. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
247001bbg	Fits Trach Sizes 00-5FR 9mm (0.354")	4.75"	1.85"	2
247002bbg	Fits Trach Sizes 6-12FR 13mm (0.512")	8"	3.66"	2

**Toothbrush Style**

Aspen Surgical offers a variety of bristle types, designed to clean all types of instruments and stains including nylon, brass, and steel. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
243001bbg	13mm (0.512"), Nylon	7.25"	1.5"	2
243002bbg	13mm (0.512"), Brass	7.25"	1.5"	2
243003bbg	13mm (0.512"), Steel	7.25"	1.5"	2

**Double-Ended Nylon Cleaning Brushes**

Our Double-Ended Cleaning Brush features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning.




Available Sterile and Non-Sterile

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity
243004bbg	Non-Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	2/Bag
243014pbx	Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	1/Pack, 12/Packs/Box

**Reusable Channel Cleaning Brushes**


Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles are available in a variety of lengths and materials to fit the range of endoscopic devices. Latex-free.



Part No.	Description	Total Length	Sheath Diameter	Quantity Per Bag
248107bbg	3-5mm	70cm	1.6mm	1
248115bbg	3-5mm	115cm	1.6mm	1
248160bbg	3-5mm	160cm	1.6mm	1
248180bbg	3-5mm	230cm	1.8mm	1
248182bbg	3-5mm — Double Ended	230cm	1.8mm	1
248185bbg	3-5mm — Teflon Coated	230cm	1.8mm	1
248190bbg	3-5mm	120cm	1.8mm	1

**Disposable Scope Channel Cleaning Brushes**

Aspen Surgical's new Disposable Scope Channel Cleaning Brushes have brush heads that range from to 2.5mm – 5mm and will fit scopes as small as 1.1mm. Our disposable scope channel cleaning brushes feature more than double the bristle length than competing low-cost brushes. They have a strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively. Single-use and latex-free.




Disposable Scope Channel Brushes

Part No.	Description	Total Length	Sheath Diameter	Quantity Per Box
248102bbg	2.5mm	230cm	1mm	50
248103bbg	3mm	230cm	1.55mm	50
248104bbg	5mm — Double Ended	230cm	1.85mm	50
248105bbg	5mm	230cm	1.85mm	50

**Hand/Nail Brushes**

Our Hand/Nail brush is used to clean dirt on hands and nails. It can also be used for general instrument cleaning. Latex-free.




Part No.	Brush Head Diameter	Quantity Per Bag
244001bbg	35mm x 93mm (1.375" x 1.75")	2

**Instrument Brush Storage**

Convenient storage with our acrylic racks helps organize your cleaning brushes and help ensure you will always have the right brush for the job. Latex-free.

Part No.	Description	Quantity
246001bbg	Brush Rack, Medium (without drip tray)	1
246003bbg	Brush Rack, Large (with drip tray)	1



In addition to Instrument Care, Aspen Surgical carries a full line O.R. Accessory, Woundcare, Patient and Staff Safety, and Ultracell® Ophthalmic and Otolaryngology products. Please call us to request one of the catalogs above, or to share your new product ideas!