

30 OTHER LAPAROSCOPIC PRODUCTS**MEDIFLEX-NUBERT TISSUE RETRIEVAL SYSTEM**

Mediflex offers the first Retrieval Bag system that features a disposable bag and reusable Forceps. The Forceps are specifically designed to insert, open, securely close and remove the Bag. The NUBERT Disposable Retrieval Bag is tough, made of tear resistant material. It is tapered to prevent bunching of tissue and it can be used repeatedly in the same procedure.

■ **REF. NO. 24003**

Mediflex-NUBERT Laparoscopic Retrieval Bag,
120ml Capacity, Box of 10, Sterile

■ **REF. NO. 91312-CFR-F**

Mediflex Modular Laparoscopic Nubert Grasper,
5mm x 33cm with On/Off Ratchet Carbon Fiber Handle

■ **REF. NO. 91312-I-F**

Replacement Insert for Mediflex-NUBERT Grasper

MEDIFLEX DISPOSABLE TISSUE RETRIEVAL SYSTEM


The Mediflex-Albert and Mediflex-Hubert bags are designed to remove larger specimens than those placed in the Nubert system. In addition the tear resistant polyurethane-coated rip-stop nylon is impervious to fluids and strong enough for use in Morcellation of tissue.

■ **REF. NO. 24001**

Mediflex-ALBERT Tissue Retrieval Bag, 900ml Capacity,
Box of 5, Sterile

■ **REF. NO. 24002**

Mediflex-HUBERT Tissue Retrieval Bag, 3000ml Capacity,
Box of 5, Sterile

**MEDIFLEX DISPOSABLE SUTURE GRASPER CLOSURE DEVICE**

Generally used for Trocar Site Closure, this versatile device can also temporarily tack hernia mesh in position, be used to assist in percutaneous suturing and to ligate abdominal wall bleeders.

A TRUE MULTI-USE INSTRUMENT

■ **REF. NO. D1001**

Sterile, 10 each per box.

