

### Clamp Covers

Clamp Covers assist in the atraumatic clamping of delicate vessels during surgery. They are available in radiopaque silicone, non-radiopaque fabric, and radiopaque fabric. We offer a variety of sizes to provide outstanding convenience and flexibility of use. These are packaged sterile or non-sterile, intended for single-use, and latex-free.



Part No.	Color	Material	Radiopaque	Length mm	Length Inches	Quantity Per Box
<b>Sterile</b>						
061011pbx	Red	Silicone	Yes	2.5 x 130	0.10 x 5	10 Pkgs of 2
061012pbx	White	Silicone	Yes	5.1 x 130	0.20 x 5	10 Pkgs of 2
061013pbx	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	10 Pkgs of 2
071001pbx	White	Fabric	No	3.2 x 130	0.125 x 5	10 Pkgs of 4
071002pbx	White	Fabric	No	9.5 x 130	0.375 x 5	10 Pkgs of 4
072001pbx	White	Fabric	Yes	3.2 x 150	0.125 x 6	10 Pkgs of 4
072002pbx	White	Fabric	Yes	9.5 x 150	0.375 x 6	10 Pkgs of 4
072004pbx	White	Fabric	Yes	6.3 x 150	0.25 x 6	10 Pkgs of 4
<b>Non-Sterile</b>						
062011bbg	Red	Silicone	Yes	2.5 x 130	0.10 x 5	20
062012bbg	White	Silicone	Yes	5.1 x 130	0.20 x 5	20
062013bbg	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	20