


Reusable Channel Cleaning Brushes


Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles are available in a variety of lengths and materials to fit the range of endoscopic devices. Latex-free.



Part No.	Description	Total Length	Sheath Diameter	Quantity Per Bag
248107bbg	3-5mm	70cm	1.6mm	1
248115bbg	3-5mm	115cm	1.6mm	1
248160bbg	3-5mm	160cm	1.6mm	1
248180bbg	3-5mm	230cm	1.8mm	1
248182bbg	3-5mm — Double Ended	230cm	1.8mm	1
248185bbg	3-5mm — Teflon Coated	230cm	1.8mm	1
248190bbg	3-5mm	120cm	1.8mm	1

Disposable Scope Channel Cleaning Brushes

Aspen Surgical's new Disposable Scope Channel Cleaning Brushes have brush heads that range from to 2.5mm – 5mm and will fit scopes as small as 1.1mm. Our disposable scope channel cleaning brushes feature more than double the bristle length than competing low-cost brushes. They have a strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively. Single-use and latex-free.




Disposable Scope Channel Brushes

Part No.	Description	Total Length	Sheath Diameter	Quantity Per Box
248102bbg	2.5mm	230cm	1mm	50
248103bbg	3mm	230cm	1.55mm	50
248104bbg	5mm — Double Ended	230cm	1.85mm	50
248105bbg	5mm	230cm	1.85mm	50

Hand/Nail Brushes

Our Hand/Nail brush is used to clean dirt on hands and nails. It can also be used for general instrument cleaning. Latex-free.




Part No.	Brush Head Diameter	Quantity Per Bag
244001bbg	35mm x 93mm (1.375" x 1.75")	2

Instrument Brush Storage

Convenient storage with our acrylic racks helps organize your cleaning brushes and help ensure you will always have the right brush for the job. Latex-free.

Part No.	Description	Quantity
246001bbg	Brush Rack, Medium (without drip tray)	1
246003bbg	Brush Rack, Large (with drip tray)	1



In addition to Instrument Care, Aspen Surgical carries a full line O.R. Accessory, Woundcare, Patient and Staff Safety, and Ultracell® Ophthalmic and Otolaryngology products. Please call us to request one of the catalogs above, or to share your new product ideas!