

Write-On ID Labels

Ideal for customization of trays, storage areas, cabinets or carts, these color-coded labels are formatted on 8½"x 11" sheets with three different label sizes on each sheet. For best results use with permanent markers only. Latex-free.



Part No.	Description
161002eea	Blue Border
161003eea	Green Border
161004eea	Red Border
161005eea	Yellow Border

Color-Coded Identification Locking Tags

These tamper evident locking tags can be customized to a variety of uses such as identifying instruments, carts, trays, containers or cabinets. Available pre-printed with "Repair", "Sharpen", or "Good" on one side and writable on the other or blank on one side and writable on the other. Sold in bags of 100. Latex-free.



Part No.	Description
251001bbg	Red, Imprint: Repair, Write-On
251002bbg	Yellow, Imprint: Sharpen, Write-On
251003bbg	Green, Imprint: Good, Write-On
251004bbg	Blue, Write-On
251005bbg	White, Write-On
251006bbg	Green, Write-On
251007bbg	Black, Write-On
251008bbg	Purple, Write-On
251009bbg	Orange, Write-On
251010bbg	Assorted Colors, Write-On

Utility Markers

Our utility markers provide a permanent method of identifying instruments, specimens, and other objects. The markers are waterproof, non-smearing, and can be used on most surfaces. Utility markers are available with or without nine black and white medication labels or twenty color medication labels and come in regular or taper tip depending on your need. Latex-free.

Regular Tip Utility Marker



Part No.	Description	Quantity Per Box
2850	Marker and 9 Label Set — Sterile	50
2850BN	Marker and 9 Label Set — Non-Sterile	100
2851	Marker — Sterile	50
2851BN	Marker — Non-Sterile	100


Taper Tip Utility Marker



Part No.	Description	Quantity Per Box
2830	Marker and 9 Label Set — Sterile	50
2830BN	Marker and 9 Label Set — Non-Sterile	100
2831	Marker — Sterile	50
2831BN	Marker — Non-Sterile	100
ML5832	Marker and 20 Medication Label Set — Sterile	50

Corner Protectors

Protect your instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays. Aspen Surgical's Corner Protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap. They can be used with steam, gas, ozone, and low temperature plasma sterilization and are reusable. Latex-free.



Part No.	Description	Quantity Per Bag
CP-80NS	Corner Protectors	80