


DISPOSABLE VERRES NEEDLES


Mediflex introduces the economic alternative to establishing Pneumoperitoneum in Laparoscopic Surgery. The Mediflex Disposable Verres Needle features an ergonomically designed non-slip handle, a sharp 14 ga. needle, a traditional spring activated blunt stylet and a clear housing to facilitate a drop test. Available in both 120 and 150mm lengths.



- REF. NO. D1002
120mm Disposable Verres Needle Sterile,
10 each per box
- REF. NO. D1003
150mm Disposable Verres Needle Sterile,
10 each per box

DISPOSABLE CHOLANGIOGRAPHY CATHETERS

Whether your surgeons introduce the Cholangiography catheter through a trocar or percutaneously ...whether they like a rigid tip or a flexible catheter, Mediflex has the right product for you.



- REF. NO. D1004
A 20" (42.8cm) flexible catheter that can be used in open or laparoscopic procedures. Complete with 14ga. percutaneous introducer. Sterile, 10 each per box.
- REF. NO. D1005
A Tapered, beaded tip 13" (27.8cm) catheter with radiolucent 5mm introducer sheath and gas seal. For laparoscopic use through a 5mm trocar/cannula. Sterile, 6 each per box.
- REF. NO. D1006
A Tapered, beaded tip 11" (23.5cm) catheter and 14ga. percutaneous introducer. Sterile, 10 each per box.

