


18 LAPAROSCOPIC INSTRUMENTS

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

Available as complete 3 piece instruments or mix and match handles, shafts and inserts for each system.



MEDIFLEX MODULAR INSTRUMENT FEATURES

	DELUXE MODULAR SYSTEM	STANDARD MODULAR SYSTEM
EASY TO CLEAN AND STERILIZE	■	■
Can be left assembled, cleaned with a flush port and sterilized as a one piece instrument	■	
Can be partially disassembled by removing the Shaft and Insert. These can be cleaned with a flush port	■	
Can be completely disassembled for Cleaning and Sterilization as a 3 piece instrument	■	■
INTERCHANGEABLE HANDLES AND SHAFTS	■	■
Ratcheted and Non-Ratcheted Handles	■	■
Dozens of inserts to choose from	■	■
Ergonomic Handle Designs	■	■
Every Pattern available in both 36cm and 45cm lengths	■	
Every Pattern available in both 33cm and 45cm lengths		■
One-handed rotation	■	■
All Handles capable of monopolar electrocautery	■	■
ALL INSERTS MADE FROM SURGICAL GRADE STAINLESS STEEL	■	■
EASY ASSEMBLY	■	■
Can be used as a one-piece instrument with no assembly	■	
Can be used as a two-piece instrument with minimal assembly	■	
INSERTS LOCKED IN PLACE WITH 2 LOCKING MECHANISMS	■	■
ECONOMICAL	■	■
Low Cost per Use	■	
Lowest Cost per Use		■

MODULAR LAPAROSCOPIC INSTRUMENT COMPONENTS

DELUXE MODULAR SYSTEM

AVAILABLE INSULATED SHAFTS WITH FLUSH PORT



■ REF. NO. 91285-D
Deluxe Insulated Shaft 5mm by 36cm



■ REF. NO. 91285-D45
Deluxe Insulated Shaft 5mm by 45cm



■ REF. NO. 91287-D
Deluxe Insulated Shaft 10mm by 36cm



■ REF. NO. 91287-D45
Deluxe Insulated Shaft 10mm by 45cm

AVAILABLE CARBON FIBER HANDLES

■ REF. NO. 91383-D-CF
Deluxe, Non-ratcheted
Ergonomic Modular Handle



■ REF. NO. 91384-D-CFR
Deluxe, On/Off Ratchet,
Ergonomic Modular Handle

STANDARD MODULAR SYSTEM

AVAILABLE INSULATED SHAFTS



■ REF. NO. 91285
Insulated Shaft 5mm by 33cm



■ REF. NO. 91285-45
Insulated Shaft 5mm by 45cm



■ REF. NO. 91287
Insulated Shaft 10mm by 33cm



■ REF. NO. 91287-45
Insulated Shaft 10mm by 45cm

AVAILABLE CARBON FIBER HANDLES

■ REF. NO. 91383-CF
Non-ratcheted Ergonomic
Modular Handle



■ REF. NO. 91384-CFR
On/Off Ratchet,
Ergonomic Modular Handle

20 LAPAROSCOPIC INSTRUMENTS

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

SCISSORS AND SHEARS

5mm Curved Metzenbaum Scissors, Flat Dissecting Blades, Double Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91467-D-I	91467-D-CF
33cm	91467-I	91467-CF
45cm	91467-D45-I	91467-D-CF-45
45cm	91467-45-I	91467-CF-45

5mm Curved Metzenbaum Scissors, Serrated Dissecting Blades, Double Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91466-D-I	91466-D-CF
33cm	91466-I	91466-CF
45cm	91466-D45-I	91466-D-CF-45
45cm	91466-45-I	91466-CF-45

5mm Straight Metzenbaum Scissors, Serrated Dissecting Blades, Double Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91465-D-I	91465-D-CF
33cm	91465-I	91465-CF
45cm	91465-D45-I	91465-D-CF-45
45cm	91465-45-I	91465-CF-45

5mm Hook Scissors, Single Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91410-D-I	91410-D-CF
33cm	91410-I	91410-CF
45cm	91410-D45-I	91410-D-CF-45
45cm	91410-45-I	91410-CF-45

5mm Micro Scissors, Straight, Single Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91344-D-I	91344-D-CF
33cm	91344-I	91344-CF
45cm	91344-D45-I	91344-D-CF-45
45cm	91344-45-I	91344-CF-45

5mm Micro Scissors, Curved, Single Action



Length	Ref. No.	
	Insert Only	with Handle & Shaft
36cm	91346-D-I	91346-D-CF
33cm	91346-I	91346-CF
45cm	91346-D45-I	91346-D-CF-45
45cm	91346-45-I	91346-CF-45

DELUXE

■ REF. NO. 91383-D-CF

Deluxe, Non-ratcheted Ergonomic Modular Handle



STANDARD

■ REF. NO. 91383-CF

Non-ratcheted Ergonomic Modular Handle



MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

DELUXE MODULAR INSTRUMENT



STANDARD MODULAR INSTRUMENT



DISSECTING FORCEPS

5mm Curved Maryland Dissector, Double Action, Cross Serrations

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91412-D-I	91412-D-CF
33cm	91412-I	91412-CF
45cm	91412-D45-I	91412-D-CF-45
45cm	91412-45-I	91412-CF-45

5mm Curved Dissector, Double Action, DeBaKey Jaws

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91317-D-I	91317-D-CF
33cm	91317-I	91317-CF
45cm	91317-D45-I	91317-D-CF-45
45cm	91317-45-I	91317-CF-45

5mm Mixer Dissector, Right Angle, Double Action, Horizontal Serrations

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91333-D-I	91333-D-CF
33cm	91333-I	91333-CF
45cm	91333-D45-I	91333-D-CF-45
45cm	91333-45-I	91333-CF-45

5mm Mixer Dissector, Right Angle, Double Action, Cross Serrations

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91376-D-I	91376-D-CF
33cm	91376-I	91376-CF
45cm	91376-D45-I	91376-D-CF-45
45cm	91376-45-I	91376-CF-45

5mm Dolphin Dissector, Straight, Double Action, Horizontal Serrations, Grasping Cup

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91415-D-I	91415-D-CF
33cm	91415-I	91415-CF
45cm	91415-D45-I	91415-D-CF-45
45cm	91415-45-I	91415-CF-45

Blue Represents **DELUXE SYSTEM**

Grey Represents **STANDARD SYSTEM**

22 LAPAROSCOPIC INSTRUMENTS

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

GRASPERS AND CLAMPS

5mm Fenestrated Grasper, 14mm Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91345-D-I	91345-D-CFR
33cm	91345-I	91345-CFR
45cm	91345-D45-I	91345-D-CFR-45
45cm	91345-45-I	91345-CFR-45

5mm Fenestrated Grasper, 15mm Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91301-D-I	91301-D-CFR
33cm	91301-I	91301-CFR
45cm	91301-D45-I	91301-D-CFR-45
45cm	91301-45-I	91301-CFR-45

5mm Fenestrated Grasper, 38mm Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91343-D-I	91343-D-CFR
33cm	91343-I	91343-CFR
45cm	91343-D45-I	91343-D-CFR-45
45cm	91343-45-I	91343-CFR-45

5mm Fenestrated Grasper, 19mm Aggressive Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91382-D-I	91382-D-CFR
33cm	91382-I	91382-CFR
45cm	91382-D45-I	91382-D-CFR-45
45cm	91382-45-I	91382-CFR-45

5mm Fenestrated Grasper, 19mm Endo-Clinch Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91323-D-I	91323-D-CFR
33cm	91323-I	91323-CFR
45cm	91323-D45-I	91323-D-CFR-45
45cm	91323-45-I	91323-CFR-45

5mm Johann Grasper, 36mm Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91311-D-I	91311-D-CFR
33cm	91311-I	91311-CFR
45cm	91311-D45-I	91311-D-CFR-45
45cm	91311-45-I	91311-CFR-45

DELUXE

■ **REF. NO. 91384-D-CFR**
Deluxe, On/Off Ratchet,
Ergonomic Carbon Fiber Handle



STANDARD

■ **REF. NO. 91384-CFR**
Standard On/Off Ratchet,
Ergonomic Carbon Fiber Handle



Blue Represents
DELUXE SYSTEM

Grey Represents
STANDARD SYSTEM

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

DELUXE MODULAR INSTRUMENT



STANDARD MODULAR INSTRUMENT



GRASPERS AND CLAMPS

5mm Johann Grasper, 56mm Double Action Jaws

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91312-D-I	91312-D-CFR
33cm	91312-I	91312-CFR
45cm	91312-D45-I	91312-D-CFR-45
45cm	91312-45-I	91312-CFR-45

5mm Johann Grasper, 76mm Double Action Jaws

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91313-D-I	91313-D-CFR
33cm	91313-I	91313-CFR
45cm	91313-D45-I	91313-D-CFR-45
45cm	91313-45-I	91313-CFR-45

5mm Maxi-Grip Fundus Grasper, 17mm Double Action Jaws

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91302-D-I	91302-D-CFR
33cm	91302-I	91302-CFR
45cm	91302-D45-I	91302-D-CFR-45
45cm	91302-45-I	91302-CFR-45

5mm Maxi-Grip Fundus Grasper, 26mm Double Action Jaws

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91377-D-I	91377-D-CFR
33cm	91377-I	91377-CFR
45cm	91377-D45-I	91377-D-CFR-45
45cm	91377-45-I	91377-CFR-45

5mm Wave Grasper, Double Action

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91375-D-I	91375-D-CFR
33cm	91375-I	91375-CFR
45cm	91375-D45-I	91375-D-CFR-45
45cm	91375-45-I	91375-CFR-45

5mm Wave Grasper with 2 x 3 Teeth, Single Action

Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91380-D-I	91380-D-CFR
33cm	91380-I	91380-CFR
45cm	91380-D45-I	91380-D-CFR-45
45cm	91380-45-I	91380-CFR-45

24 LAPAROSCOPIC INSTRUMENTS

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

GRASPERS AND CLAMPS

5mm Standard Grasper, Full Serrations, Double Action Jaws



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91354-D-I	91354-D-CFR
33cm	91354-I	91354-CFR
45cm	91354-D45-I	91354-D-CFR-45
45cm	91354-45-I	91354-CFR-45

5mm Allis Tissue Forceps, Rounded Tip, Double Action Jaws



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91337-D-I	91337-D-CFR
33cm	91337-I	91337-CFR
45cm	91337-D45-I	91337-D-CFR-45
45cm	91337-45-I	91337-CFR-45

5mm Babcock Tissue Forceps, Double Action DeBakey Jaws



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91542-D-I	91542-D-CFR
33cm	91542-I	91542-CFR
45cm	91542-D45-I	91542-D-CFR-45
45cm	91542-45-I	91542-CFR-45

10mm Babcock Tissue Forceps, Double Action DeBakey Jaws



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91545-D-I	91545-D-CFR
33cm	91545-I	91545-CFR
45cm	91545-D45-I	91545-D-CFR-45
45cm	91545-45-I	91545-CFR-45

10mm Atraumatic Spring Jaw Babcock Tissue Forceps, Double Action



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91544-D-I	91544-D-CFR
33cm	91544-I	91544-CFR
45cm	91544-D45-I	91544-D-CFR-45
45cm	91544-45-I	91544-CFR-45

5mm DeBakey Grasper, 60mm Double Action Jaws



Length	Ref. No.	Ref. No.
	Insert Only	with Handle & Shaft
36cm	91565-D-I	91565-D-CFR
33cm	91565-I	91565-CFR
45cm	91565-D45-I	91565-D-CFR-45
45cm	91565-45-I	91565-CFR-45

DELUXE

■ **REF. NO. 91384-D-CFR**
Deluxe, On/Off Ratchet,
Ergonomic Carbon Fiber Handle



STANDARD

■ **REF. NO. 91384-CFR**
Standard On/Off Ratchet,
Ergonomic Carbon Fiber Handle



Blue Represents
DELUXE SYSTEM

Grey Represents
STANDARD SYSTEM

MODULAR LAPAROSCOPIC INSTRUMENT SYSTEMS

DELUXE MODULAR INSTRUMENT



STANDARD MODULAR INSTRUMENT



TENACULA, CLAWS AND SPOON FORCEPS

5mm Tenaculum, Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91537-D-I	91537-D-CFR
33cm	91537-I	91537-CFR
45cm	91537-D45-I	91537-D-CFR-45
45cm	91537-45-I	91537-CFR-45

10mm Tenaculum, Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91536-D-I	91536-D-CFR
33cm	91536-I	91536-CFR
45cm	91536-D45-I	91536-D-CFR-45
45cm	91536-45-I	91536-CFR-45

5mm Claw Forceps, Single Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91606-D-I	91606-D-CFR
33cm	91606-I	91606-CFR
45cm	91606-D45-I	91606-D-CFR-45
45cm	91606-45-I	91606-CFR-45

10mm Claw Forceps, Single Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91605-D-I	91605-D-CFR
33cm	91605-I	91605-CFR
45cm	91605-D45-I	91605-D-CFR-45
45cm	91605-45-I	91605-CFR-45

5mm Spoon/Stone Retrieval Forceps, Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91485-D-I	91485-D-CFR
33cm	91485-I	91485-CFR
45cm	91485-D45-I	91485-D-CFR-45
45cm	91485-45-I	91485-CFR-45

10mm Spoon/Stone Retrieval Forceps, Double Action Jaws



Length	Ref. No. Insert Only	Ref. No. with Handle & Shaft
36cm	91660-D-I	91660-D-CFR
33cm	91660-I	91660-CFR
45cm	91660-D45-I	91660-D-CFR-45
45cm	91660-45-I	91660-CFR-45

26 LAPAROSCOPIC INSTRUMENTS**ONE-PIECE, FLUSHABLE LAPAROSCOPIC INSTRUMENTS****10mm ANVIL GRASPER WITH FLUSH PORT**

■ REF. NO. 91599-A

10mm x 33cm Anvil Grasper
with Flush Port.

■ REF. NO. 91599-A-45

10mm x 45cm Anvil Grasper
with Flush Port.

**ARTICULATING FAN RETRACTORS WITH FLUSH PORT**

■ REF. NO. 91671

5mm Articulating Fan Retractor,
3 blades with Flush Port.



■ REF. NO. 91672

10mm Articulating Fan Retractor,
5 blades with Flush Port.

**ARTICULATING ESOPHAGEAL DISSECTOR/RETRACTOR WITH FLUSH PORT**

■ REF. NO. 91667

5mm x 33cm Articulating Esophageal Dissector and Retractor

■ REF. NO. 91667-45


5mm x 45cm Articulating Esophageal Dissector and Retractor

■ REF. NO. 91668

10mm x 33cm Articulating Esophageal Dissector and Retractor

■ REF. NO. 91668-45

10mm x 45cm Articulating Esophageal Dissector and Retractor

**LAPAROSCOPIC NEEDLE HOLDERS WITH FLUSH PORT**

■ REF. NO. 91592-AX


Needle Holder, One-piece with Flush Port,
5mm x 33cm, Axial Handle, Straight Jaws

■ REF. NO. 91593-AX

Needle Holder, One-piece with Flush Port,
5mm x 33cm, Axial Handle, Jaws Curve Right

■ REF. NO. 91594-AX

Needle Holder, One-piece with Flush Port,
5mm x 33cm, Axial Handle, Jaws Curve Left



OTHER LAPAROSCOPIC INSTRUMENTATION
AUTOCLAVABLE LAPAROSCOPES



■ **REF. NO. 91120**
5mm x 30.5cm
0° View

■ **REF. NO. 91125**
5mm x 30.5cm
30° View

■ **REF. NO. 91130**
10mm x 31.5cm
0° View

■ **REF. NO. 91135**
10mm x 31.5cm
30° View

■ **REF. NO. 91136**
10mm x 31.5cm
45° View




REUSABLE KITE STYLE FAN RETRACTORS

■ **REF. NO. 91670**
5mm x 33cm Straight Kite Retractor

■ **REF. NO. 91675**
5mm x 33cm Curved Kite Retractor

REUSABLE MONOPOLAR ELECTRODES



SPATULA

L-HOOK

J-HOOK

■ **REF. NO. 91170**
5mm x 33cm Spatula Tip Electrode

■ **REF. NO. 91175**
5mm x 33cm L-Hook Tip Electrode

■ **REF. NO. 91185**
5mm x 33cm J-Hook Tip Electrode

■ **REF. NO. 91170-45**
5mm x 45cm Spatula Tip Electrode

■ **REF. NO. 91175-45**
5mm x 45cm L-Hook Tip Electrode

■ **REF. NO. 91185-45**
5mm x 45cm J-Hook Tip Electrode

REUSABLE "S" RETRACTORS FOR HASSON OPEN LAPAROSCOPY




■ **REF. NO. 69080**
"S" Retractors 1/4" Wide (Sold as Pair)

■ **REF. NO. 69081**
"S" Retractors 3/8" Wide (Sold as Pair)

■ **REF. NO. 69082**
"S" Retractors 1/2" Wide (Sold as Pair)

REUSABLE MULTI-PURPOSE KNOT PUSHER/PALPATION PROBE



ONE END IS A PALPATION PROBE

Detail of Knot Pusher

■ **REF. NO. 94435**
4mm x 36cm Multi-Purpose Knot Pusher

■ **REF. NO. 94436**
4mm x 45cm Multi-Purpose Knot Pusher

THE OTHER END IS A KNOT PUSHER

28 LAPAROSCOPIC INSTRUMENT HOLDING/POSITIONING**THE MEDIFLEX-DRAUGHN LAP-RACK
INSTRUMENT HOLDING AND POSITIONING SYSTEM****QUICK AND EASY ACCESS TO LAPAROSCOPIC INSTRUMENTS
WHEN AND WHERE YOU NEED THEM!**

The Lap Rack Instrument Holding and Positioning System is designed to hold up to eight laparoscopic instruments of any size including clip applicators, graspers, scissors, and stapling devices. It can also enhance cord management by holding cautery probes, suction irrigators and harmonic scalpels. Developed by David Draughn, MD of Gastonia, NC.

**Cantilever-Style Rack Attaches Directly
to OR Table Rail**

The lap rack system moves with the surgical table. Allows positioning exactly where the surgeon needs the instruments.

**Eliminates Awkward Hand-to-Hand
Instrument Exchanges**

The surgeon chooses the correct instrument at his or her convenience. The scrub tech no longer has to "guess" what's needed next.

**Keeps Instruments In the Sterile Field—
Not On the Floor**

Virtually eliminates disruption of procedure to replace instruments.

**Provides Immediate Access to Instruments
Without Losing Sight of Monitor**

Surgeon maintains focus on the procedure.

Can Save Time and Labor

Procedures with fewer interruptions save valuable OR time.

**ALSO FOR USE WITH THE LAP RACK SYSTEM:
FLEX ARM™ COMPANION ARM**

Useful for holding laparoscopic devices such as scopes, graspers and retractors. 69303-AQC (shown at right) attaches to the horizontal or vertical portion of the Mount Post.

Please see page 8 for additional information.



■ REF. NO. 69303-AQC

■ REF. NO. 69074

The Mediflex Draughn Lap Rack System Complete with:

Qty.	Ref. No.	Description/Size
1 ea.	Ref. No. 69097-FR	Lap Rack and Lap Rack Mount
1 ea.	Ref. No. 69098	Long Reach Mount Post/ Swivel Ball Assembly
1 ea.	Ref. No. 69099	OR Rail Clamp

CUTTING EDGE™ REPOSABLE SCISSOR AND DISSECTOR TIPS

The Cutting Edge Product Line from Mediflex is a Laparoscopic System consisting of Reusable Laparoscopic Instrument Handles and Disposable Scissor and Dissector tips. It provides the economy of a reusable handle with the guarantee of a new tip for each procedure. Photos below show 2 views of each tip: opened, closed from above and closed from the side.



■ **REF. NO. 91710-10**
Curved 3/4" (1.9cm)
Metzenbaum Scissor Tips,
Box of 10, Sterile



■ **REF. NO. 91740-10**
Hook Scissor Tips, Box
of 10, Sterile



■ **REF. NO. 91720-10**
Straight 5/8" (1.6cm)
Metzenbaum Scissor Tips,
Box of 10, Sterile



■ **REF. NO. 91750-10**
Curved 3/4" (1.9cm)
Maryland Dissector Tips,
Box of 10, Sterile




■ **REF. NO. 91730-10**
Mini Curved 1/2"
(1.25cm) Metzenbaum
Scissor Tips, Box of 10,
Sterile



■ **REF. NO. 91760-10**
Baby Curved 3/8" (.8mm)
Metzenbaum Scissor Tips,
Box of 10, Sterile

AVAILABLE REUSABLE HANDLES CUTTING EDGE SYSTEM



■ **REF. NO. 91383-45TI**
45cm x 5mm Handle, Complete with Flush
Port, Monopolar Post and Rotation

■ **REF. NO. 91383TI**
33cm x 5mm Handle, Complete with Flush
Port, Monopolar Post and Rotation

■ **REF. NO. 91383-26TI**
26cm x 5mm Handle, Complete with Flush
Port, Monopolar Post and Rotation

30 OTHER LAPAROSCOPIC PRODUCTS**MEDIFLEX-NUBERT TISSUE RETRIEVAL SYSTEM**

Mediflex offers the first Retrieval Bag system that features a disposable bag and reusable Forceps. The Forceps are specifically designed to insert, open, securely close and remove the Bag. The NUBERT Disposable Retrieval Bag is tough, made of tear resistant material. It is tapered to prevent bunching of tissue and it can be used repeatedly in the same procedure.

■ **REF. NO. 24003**

Mediflex-NUBERT Laparoscopic Retrieval Bag,
120ml Capacity, Box of 10, Sterile

■ **REF. NO. 91312-CFR-F**

Mediflex Modular Laparoscopic Nubert Grasper,
5mm x 33cm with On/Off Ratchet Carbon Fiber Handle

■ **REF. NO. 91312-I-F**

Replacement Insert for Mediflex-NUBERT Grasper

MEDIFLEX DISPOSABLE TISSUE RETRIEVAL SYSTEM


The Mediflex-Albert and Mediflex-Hubert bags are designed to remove larger specimens than those placed in the Nubert system. In addition the tear resistant polyurethane-coated rip-stop nylon is impervious to fluids and strong enough for use in Morcellation of tissue.

■ **REF. NO. 24001**

Mediflex-ALBERT Tissue Retrieval Bag, 900ml Capacity,
Box of 5, Sterile

■ **REF. NO. 24002**

Mediflex-HUBERT Tissue Retrieval Bag, 3000ml Capacity,
Box of 5, Sterile


**MEDIFLEX DISPOSABLE SUTURE GRASPER CLOSURE DEVICE**

Generally used for Trocar Site Closure, this versatile device can also temporarily tack hernia mesh in position, be used to assist in percutaneous suturing and to ligate abdominal wall bleeders.

A TRUE MULTI-USE INSTRUMENT


■ **REF. NO. D1001**

Sterile, 10 each per box.



DISPOSABLE VERRES NEEDLES

Mediflex introduces the economic alternative to establishing Pneumoperitoneum in Laparoscopic Surgery. The Mediflex Disposable Verres Needle features an ergonomically designed non-slip handle, a sharp 14 ga. needle, a traditional spring activated blunt stylet and a clear housing to facilitate a drop test. Available in both 120 and 150mm lengths.



- **REF. NO. D1002**
120mm Disposable Verres Needle Sterile,
10 each per box
- **REF. NO. D1003**
150mm Disposable Verres Needle Sterile,
10 each per box

DISPOSABLE CHOLANGIOGRAPHY CATHETERS

Whether your surgeons introduce the Cholangiography catheter through a trocar or percutaneously ...whether they like a rigid tip or a flexible catheter, Mediflex has the right product for you.



- **REF. NO. D1004**
A 20" (42.8cm) flexible catheter that can be used in open or laparoscopic procedures. Complete with 14ga. percutaneous introducer. Sterile, 10 each per box.
- **REF. NO. D1005**
A Tapered, beaded tip 13" (27.8cm) catheter with radiolucent 5mm introducer sheath and gas seal. For laparoscopic use through a 5mm trocar/cannula. Sterile, 6 each per box.
- **REF. NO. D1006**
A Tapered, beaded tip 11" (23.5cm) catheter and 14ga. percutaneous introducer. Sterile, 10 each per box.

